

**900ES**  
**Liter**  
Aluminum

**ecosmart™ FREE Return Shipping Instructions:**

1. Before returning cylinder, twist valve open. This will release any excess pressure.
2. PEEL OFF label wrap by pulling on the perforated strip. (Empty sticker is shown)
3. Generate your FREE label by using the ecosmart returns button at [www.gascogas.com](http://www.gascogas.com)
4. Prepare a box with adequate packing material and ship FREE to GASCO!

**Material**  
Aluminum Alloy

**Water Capacity**  
4.78 Liters

**Pressure**  
2000 psi

**Weight**  
15.8 lbs

**Valve Inlet**  
CGA 330, CGA 350,  
CGA 580, CGA 590,  
CGA 660

**Gas Type**  
Reactive  
Non-reactive

**Dimensions**  
19.25" x 5.25"



**GASCO**

— Calibration Gas —

Setting Standards. Exceeding Expectations.

[www.gascogas.com](http://www.gascogas.com) | 1-800-910-0051



PRECISION GAS MIXTURE	
COMPONENTS	
CONCENTRATION (MOLES)	
CARBON MONOXIDE (CO)	100 PPM
HYDROGEN SULFIDE (H2S)	25 PPM
PENTANE (C5H12)	0.35% VOL (25% LEL)
OXYGEN (O2)	19% VOL
NITROGEN (N2)	BALANCE

**ecosmart™ REFILLABLE CYLINDER**  
**DO NOT DISCARD**

**INSTRUCTIONS ON RETURNING ecosmart CYLINDERS**  
1. Install ecoM™ plug to ensure cylinder is empty. (COI) completed & reusable can also be used.  
2. PEEL OFF label wrap by pulling on perforated strip (empty sticker is shown).  
3. For free return shipping label go to [www.gascogas.com](http://www.gascogas.com) and click on the "ecosmart returns" button.  
4. Complete form. Labels will be generated and available to print through links at the top of the page.  
5. Ship via FEDEX with adequate and reasonable packing material.

Unlike traditional calibration gas cylinders, ecosmart was made to be returned. Simply remove the label wrap to reveal return shipping instructions. It's that easy.

**RETURN**